OWANDY-105

COMMUNICATION MADE EASY

WIRELESS Instructions manual Rev01

Intraoral Scanner

Instructions, Gebrauchsanweisung,
Instrucciones, Instructions, Istruzioni

Contents


| | |
|-----------------------------------|-----------|
| Labels and Symbols | 5 |
| EN/English..... | 9 |
| Notice | 9 |
| Description..... | 10 |
| Setting Up the Scanner | 17 |
| Getting Started | 22 |
| Safety Guide..... | 26 |
| Troubleshooting | 34 |
| Regulatory Information..... | 36 |
| Technical Specifications..... | 42 |
| DE/Deutsch..... | 45 |
| Hinweis | 45 |
| Beschreibung..... | 46 |
| Einrichtung | 54 |
| Erste Schritte | 59 |
| Sicherheitsleitfaden | 63 |
| Fehlersuche | 72 |
| Regulatorische Informationen..... | 74 |
| Technische Daten | 81 |
| ES/ Español | 83 |
| Aviso | 83 |
| Descripción | 84 |
| Configuración del Escáner | 91 |
| Primeros pasos | 96 |
| Guía de seguridad..... | 100 |


| | |
|-----------------------------------|------------|
| Solución de problemas | 109 |
| Información reglamentaria | 111 |
| Especificaciones técnicas..... | 118 |
| FR/Français | 120 |
| Avis | 120 |
| Description..... | 121 |
| Configuration du Scanner | 128 |
| Pour commencer | 133 |
| Guide de sécurité..... | 137 |
| Dépannage..... | 146 |
| Informations réglementaires | 148 |
| Spécifications techniques | 155 |
| IT/Italiano | 158 |
| Avviso | 158 |
| Descrizione | 159 |
| Impostazione del Scanner | 166 |
| Come iniziare | 171 |
| Guida alla sicurezza..... | 175 |
| Risoluzione dei problemi | 184 |
| Informazioni normative | 186 |
| Specifiche tecniche | 193 |
| Contact Information..... | 195 |

Labels and Symbols

| | |
|---|---|
|  | <p>EN Type BF applied part</p> <p>DE Anwendungsteil vom Typ BF</p> <p>ES Componente aplicado de tipo BF</p> <p>FR Partie appliquée de type BF</p> <p>IT Parte applicata di tipo BF</p> |
|  | <p>EN Separate collection for waste electric and electronic equipment (WEEE) is required.</p> <p>DE Die getrennte Sammlung von Elektro- und Elektronikgeräten (WEEE) ist vorgeschrieben.</p> <p>ES Es necesario recoger por separado los residuos de aparatos eléctricos y electrónicos (RAEE).</p> <p>FR La collecte sélective des déchets d'équipements électriques et électroniques (DEEE) est obligatoire.</p> <p>IT È richiesta la raccolta differenziata dei rifiuti di apparecchiature elettriche ed elettroniche (RAEE).</p> |
|  | <p>EN Manufacturer</p> <p>DE Hersteller</p> <p>ES Fabricante</p> <p>FR Fabricant</p> <p>IT Fabbricante</p> |
|  | <p>EN Date of Manufacture</p> <p>DE Herstellungsdatum</p> <p>ES Fecha de manufactura</p> <p>FR Date de fabrication</p> <p>IT Data di Produzione</p> |
|  | <p>EN Importer</p> <p>DE Importeur</p> <p>ES Importador</p> <p>FR Importateur</p> <p>IT Importatore</p> |
|  | <p>EN Refer to instruction manual/booklet</p> <p>DE Siehe Gebrauchsinformation</p> <p>ES Consultar el manual de instrucciones/manual</p> <p>FR Se reporter au manuel/notice d'utilisation</p> <p>IT Consultare le istruzioni del Manuale d'Uso</p> |

| | |
|---|--|
|  | <p>EN Direct current</p> <p>DE Gleichstrom</p> <p>ES Corriente continua</p> <p>FR Courant continu</p> <p>IT Corrente continua</p> |
|  | <p>EN Sterilizable in a steam sterilizer (autoclave) at 134°C</p> <p>DE Sterilisierbar in einem Dampfsterilisator (Autoklav) bei 134°C</p> <p>ES Esterilizable en esterizador de vapor (autoclave) a 134 °C</p> <p>FR Stérilisable dans un stérilisateur à vapeur (autoclave) à 134°C</p> <p>IT Sterilizzabile in uno sterilizzatore a vapore (autoclave) a 134°C</p> |
|  | <p>EN Model Number</p> <p>DE Modellnummer</p> <p>ES Número de modelo</p> <p>FR Numéro de modèle</p> <p>IT Numero di modello</p> |
|  | <p>EN Serial number</p> <p>DE Seriennummer</p> <p>ES Número de serie</p> <p>FR Numéro de série</p> <p>IT Numero di serie</p> |
|  | <p>EN Rechargeable battery</p> <p>DE Wiederaufladbare Batterie</p> <p>ES Batería recargable</p> <p>FR Batterie rechargeable</p> <p>IT Batteria ricaricabile</p> |
|  | <p>EN CE marking applicable for European Union</p> <p>DE CE-Kennzeichnung gilt für die Europäische Union</p> <p>ES Marca CE aplicable para la Unión Europea</p> <p>FR Marquage CE applicable à l'Union européenne</p> <p>IT Marchio CE applicabile all'Unione Europea</p> |
|  | <p>EN This device complies with Part 15 of the FCC Rules</p> <p>DE Dieses Gerät entspricht Teil 15 der FCC-Vorschriften</p> <p>ES Este dispositivo cumple con la parte 15 de las normas de la FCC</p> |

| | |
|---|--|
| | <p>FR Cet appareil est conforme à la partie 15 des règles de la FCC.</p> <p>IT Questo dispositivo è conforme alla Parte 15 delle norme FCC.</p> |
|  | <p>EN Authorized Representative in the European Community</p> <p>DE Bevollmächtigter in der Europäischen Gemeinschaft</p> <p>ES Representante autorizado en la Comunidad Europea</p> <p>FR Mandataire Européen</p> <p>IT Rappresentante autorizzato per la Comunità Europea</p> |
|  | <p>EN Medical Device</p> <p>DE Medizinisches Gerät</p> <p>ES Productos sanitarios</p> <p>FR Dispositif médical</p> <p>IT Dispositivo medico</p> |
|  | <p>EN Atmospheric pressure limitation</p> <p>DE Atmosphärische Druckbegrenzung</p> <p>ES Limitación de la presión atmosférica</p> <p>FR Limitation de la pression atmosphérique</p> <p>IT Limitazione della pressione atmosferica</p> |
|  | <p>EN Fragile, handle with care</p> <p>DE Zerbrechlich, mit Vorsicht zu behandeln</p> <p>ES Frágil, manipular con cuidado</p> <p>FR Fragile, à manipuler avec précaution</p> <p>IT Fragile, da maneggiare con cura</p> |
|  | <p>EN Humidity limitation</p> <p>DE Begrenzung der Luftfeuchtigkeit</p> <p>ES Limitación de la humedad</p> <p>FR Limitation de l'humidité</p> <p>IT Limitazione dell'umidità</p> |
|  | <p>EN Keep away from rain</p> <p>DE Vom Regen fernhalten</p> <p>ES Mantener alejado de la lluvia</p> <p>FR Tenir à l'écart de la pluie</p> <p>IT Tenere lontano dalla pioggia</p> |

| | |
|---|---|
|  | <p>EN Stacking limit by number</p> <p>DE Stapelungsgrenze nach Anzahl</p> <p>ES Límite de apilamiento por número</p> <p>FR Limite d'empilage par nombre</p> <p>IT Limite di impilamento per numero</p> |
|  | <p>EN Temperature limits</p> <p>DE Temperaturgrenzwerte</p> <p>ES Límites de temperatura</p> <p>FR Limites de température</p> <p>IT Limiti di temperatura</p> |
|  | <p>EN This way up</p> <p>DE Dieser Weg nach oben</p> <p>ES Por aquí arriba</p> <p>FR Par ici</p> <p>IT Da qui in su</p> |

EN/English

Notice

This document includes information on the safety instructions, regulatory information, Technical Specifications and installation and operation instructions of the device. We recommend that you thoroughly familiarize yourself with this document to make the most effective use of your system.

The information in this document is subject to change. Neither the provider nor any of its subsidiaries shall be liable for errors contained herein or for incidental damages in conjunction with the furnishing, performance, or use of this material. No part of this publication may be reproduced without the permission of the provider.

Conventions

The following special messages emphasize information or indicate potential risks to personnel or equipment.

| | |
|---|--|
|  | WARNING: Avoid injury to yourself or others by following the safety instructions precisely. |
|  | CAUTION: Alerts you to a condition that might cause serious damage or cause problems. |
| | NOTE: Provides extra information and hints. |

WARNING: We recommend that you consult this document before using the device.

All trademarks and registered trademarks are the property of their respective holders.

The device is intended for professional use only.

U.S. Federal law restricts this device to sale by or on the order of a dentist.

If any serious incident occurs in relation to the device, the user must report it to the provider and to the competent authority of its Member State in the European Union.

Description

Intended used

The **Intraoral Scanner** (herein after referred as '**scanner**') is a digital optical scanning device used to obtain digital impressions of hard and soft tissues such as teeth, gums, and mucous membrane by means of oral scanning, for oral restoration and orthodontic treatment of malocclusion.

The scanner could be used for both adults and children in clinical practice.

Contraindications

None


Clinical Benefits and Performance Characteristics


The scanner benefits dental practices by enabling practitioners to acquire digital impressions with the quality and accuracy required for digital CAD/CAM dental applications. The actual performance of the device is dependent on the user's training and operating execution. The user is solely responsible for the accuracy, completeness, and adequacy of the acquired data.

The scanner is designed to acquire 3D models in the followings:

- Upper jaw
- Lower jaw
- Buccal bite registration

Nomenclature

| |
|---|
| [1] Handpiece |
| [2] Scanner Tip |
| [3] Scan Button |
| <ul style="list-style-type: none"> ● Press once to start scanning, press again to stop scanning |
| [4] Mode Indicator |
|  Upper jaw scan mode |
|  Lower jaw scan mode |
|  Buccal bite registration mode |
| [5] Mode Switch Button |
| <ul style="list-style-type: none"> ● Press this button to switch between different modes |

| |
|---|
| [6] Status Indicator |
| <ul style="list-style-type: none">● Blinking slowly: starting up● Blinking rapidly: preparing network● Breathing: ready to be connected to the software● On: Connected to the software and selected as active scanner |
| [7] Battery Indicator |
| <ul style="list-style-type: none">● Solid green: battery level higher than 50%● Solid yellow: battery level between 20% and 50%● Blinking yellow: battery level lower than 20%● Breathing green: charging and battery level higher than 50%● Breathing yellow: charging and battery level lower than 50%● Solid green: fully charged |
| [8] Lens Window |
| [9] Charging Holder |
| [10] USB Cable (1m length) |

Computer System Requirements


For the computer system requirements, refer to the 'Technical Specifications' section


WARNING: It is MANDATORY to check that your computer system configuration is compatible with the computer system requirements for the accompanying software.

Charging the Battery

A fully charged battery can support up to 2 hours of scan time. When the Battery Indicator flashes yellow, it means the battery is low and needs to be charged in time. To charge, the scanner needs to be inserted into the Charging Holder.


The Charging Holder must be connected to a USB Type C power source to be functional.


You can also continue to use the scanner by replacing the battery, follow these steps to replace the battery.

- Turn the battery cover with the small round piece in the wall mounting kit or a coin to unlock.


- Take out the low battery and insert a charged one.


WARNING: Please use the battery specified by the provider (Li-18650-3.6V 3500mAh-PCM-Cap)

- Close the battery cover by turning it to lock.

- The charging holder holds a spare battery. Open the battery cover by turning it to unlock and you will be able to take the battery out.


NOTE: The spare batteries inside the charging holder are not charged when the scanner is in the charging holder (the spare batteries inside the charging holder are only charged when the scanner is removed from the charging holder).

NOTE: Please fully charge the battery before use.


Wireless Connection

The scanner connects to the software through a wireless connection. Before using the scanner, a compatible wireless adapter needs to be connected to the computer that runs the software.

NOTE: Expand the wireless adapter to maximize the network signal strength.


After the compatible wireless adapter is connected, Select the scanner you want to connect from the Scanner list window at the bottom right of the software and click the connect button.

| | |
|---|---|
|  |  |
| <p>Indicates a wireless scanner</p> | <p>Indicates a wired scanner</p> |

Charging Holder Overview

Place the scanner in the Charging Holder when you are not using it.


NOTE: The scanner will automatically shut down if it is idle for more than 3 minutes. If you need to use it again, press any button on the scanner, or take it out of the Charging Holder.

Setting Up the Scanner

Setting Up the Scanner

To set up the scanner, follow these steps:

- Visit the website to download the software installation package according to the product model (please refer to the download instruction page in the suitcase for website address).
- Double-click the software installation file.
- Choose a language from the drop-down list and click Ok to install.
- Follow the instructions on the screen to complete the installation.
- Firmly slide the Scanner Tip onto the end of the scanner.

- Short press any button to start the scanner.
- After the scanner is powered on, open the software and connect the scanner through the wireless network.


Using the Scanner Charging Holder

The Charging Holder can be used standing on a flat plane or mounted on the wall.


Installing the Desktop Charging Holder

To place the Charging Holder on a desktop, follow these steps:


- Select a clean surface area.
- Insert one end of the charging cable into the bottom of the charger.


- Insert the charger with the connected charging cable into the base.

- Insert the charging cable into the holder bracket.

- Pass the other end of the cable through the hole reserved on the base.


- Insert the scanner into the Charging Holder.

Installing the Wall Mount Charging Holder

To mount the Charging Holder on the wall, follow these steps:


- Select an area that you can access easily.


- Insert the screws through the holes in the holder bracket to affix it to a solid surface.
- Insert one end of the charging cable into the bottom of the charger.

- Insert the Charging Holder into the holder bracket.

- Insert the scanner into the Charging Holder.

 **CAUTION:** If the Charging Holder is not properly installed, there is a risk that the Charging Holder can fall off the wall, resulting in damage to the scanner.

Getting Started

Accessing the software user interface

To access the software user interface, follow these steps:

- Double-click the software icon on the desktop.
- The software will automatically pop up the login window.
- If you don't have an account, click **Create an account** to register the organization, and complete the email verification.
- Type your account information in the username and password fields, click the **Log in** button.
- The **Cases** page will be shown up.
- Click the **New Case** button, then input patient name and click the **OK** button.
- Click the **Scan** button.
- If the scanner is not activated, the device activation dialog will be displayed. Follow the instructions on the screen to complete device activation.
- Click the Option menu button  and select **Preferences**.
- Customize the configuration options.
- You can now start acquiring 3D models.

Audio Prompts

Scanning sound

When the **Enable scanning sound** option is enabled, your computer will play a continuous sound when you are successfully scanning. If the sound stops, it means the scan has stopped. If you need to continue, please return to the previous scanning area until the scanner resumes scanning and your computer plays a continuous sound. When a bite registration relationship is successfully scanned, your computer will also play a short sound.

Warning sound

When the **Enable warning Sound** option is enabled, if the cumulative scanning time of the current case exceeds the recommended threshold and your computer may not be able to maintain the peak scanning performance, your computer will play a short warning sound.

NOTE: If your computer does not have speakers, these options will not take effect.

Preparing the Teeth

- If there is a preparation area, retract the gingiva by gingival restriction cords. And extract the cords just before scanning the preparation.
- Before starting the scan, dry the teeth thoroughly.
- During the scan, re-dry the teeth moderately.

Preparing the Scanner


The Scanner Tip attached to the Handpiece provides sanitary shield for patients. Always disinfect the body of the scanner and perform sterilization on the Scanner Tip after each use.

CAUTION: Scanner Tips are shipped non-sterile. You must sterilize them before the first use.


CAUTION: Avoid any liquid from leaking into the air outlet near the Scanner Tip mount or the air inlet at the rear of the Handpiece, otherwise the Handpiece may be damaged.

To prepare the scanner, follow these steps:

- Make sure the Lens Window at the base of the scanner is clean by wiping it with a moist, lint-free cloth or lens tissue.
- Slide the Scanner Tip onto the scanner as shown below.


Starting Scanning

To start scanning, place the Scanner Tip on the surface of the tooth to stabilize the scanner and press the Scan button. Wait until a 3D model appears on the 3D model display screen, and then slowly move it along the arch at 0-5mm from the teeth.

Scanning Approach

The recommended scanning method is to start with a molar, since it has greater details for easier identification. Change the scanning angle to less than 60 degrees during scanning to allow the surfaces to overlap, if the overlap is too small, the alignment may be lost.


Scanning Protocol

The recommended scanning protocol consisting of 3 sweeps: occlusal, lingual and buccal to ensure good data coverage of all surfaces.

It is recommended to start the first sweep from the occlusal surface, and you should start with the first molar.

The second sweep can scan both the lingual and buccal sides, and the third scan covers the opposite side of the second sweep.

Safety Guide

Warnings and Safety Instructions

DANGER OF ELECTRIC SHOCK

This is an electrical unit. Do NOT expose it to water spray. Such action can cause an electric shock or a malfunction of the unit.

CAUTION: All known residual risks, contraindications, or undesirable side effects are listed in this document. If any serious incident occurs in relation to the device, you must report it to the provider and to the competent authority of your Member State in the European Union.

WARNING

Scanner

- You **MUST** read and understand this safety information before using the scanner.
- This scanner shall only be used inside hospitals and other professional healthcare facilities and **MUST NOT** be used near high frequency surgical equipment and the RF shielded room of an ME System for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.
- Before using the scanner, check the outer surfaces of the unit and any accessories to ensure there are no rough surfaces, sharp edges, or protrusions which may cause a safety hazard.
- You are responsible for the operation and maintenance of the scanner. You **MUST** have training to use the scanner.
- **DO NOT** place objects within the field of operation of the unit.
- When the unit is not in use, ensure that the scanner is turned **OFF**.
- **DO NOT** use the scanner in conjunction with oxygen-rich environments. This unit is not intended for use with flammable anesthetics or flammable agents.
- **DO NOT** pull or twist the cable.
- **DO NOT** drop the scanner or the accessories.
- **DO NOT** sterilize the scanner.

- DO NOT expose the scanner to water spray or submerge it in water or disinfectant.
- DO NOT expose the scanner to high vibrations.
- DO NOT expose the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- DO NOT stare at the LED emission window.
- When the Scanner Tip is removed, install the front protective cover to protect the Lens Window.
- DO NOT remove the cover of any scanner components. The scanner contains no user-serviceable parts. For any repairs, contact a qualified the provider service technician.
- DO NOT replace the cables provided with the scanner with other cables. Doing so may damage the scanner and adversely affect the safety protection and EMC performance of the scanner.
- Any other equipment not complying with IEC-60601 shall be kept at least 1.5 meters away from the patient.
- If the equipment is faulty, turn it OFF, display an “Out of Service” notice, and contact a qualified the provider service technician.
- Using components, accessories, cables, and spare parts other than those specified or provided by the manufacturer of this equipment may impair the safety protection of the scanner and may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- No modification of this equipment is allowed.
- Additional multiple outlet strips or extension cords should not be connected to the system.
- The maximum temperature of the applied part may reach to 43 °C; to avoid overheating, do not use it for extended periods.
- To isolate the holder from power supply, unplug the USB connector from the USB port.
- DO NOT maintain or service this equipment while it is in use with the patient.
- Connection of the PEMS (Programmable Electrical Medical System) to an IT NETWORK that includes other equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.
- Patients with oral mucosal disease, mental illness, severe respiratory disease, asthma, Parkinson's disease, hyperactivity disease are forbidden.
- Patients with moderate or severe opening limitation should use it with caution.

Computer

- Do NOT place any equipment which does not comply with IEC 60601-1 in the immediate vicinity of the patient.

Leave at least 1.5 meters distance between the patient and the equipment.

- The scanner is only intended to be connected to a computer that is at least IEC 60950 / IEC 62368, or equivalent standards certified. Connecting the scanner to other equipment may be hazardous.
- See the installation guide for your computer for information about the data processing system, computer, and screen. Leave a sufficient amount of clear space around the computer to ensure that it is properly ventilated.
- Position the screen to avoid light reflections from internal or external lighting for maximum image quality and visual comfort.

Battery

- Do not dismantle, open, or shred secondary cells or batteries.
- Do not expose cells or batteries to heat or fire. Avoid storage in direct sunlight.
- Do not short-circuit a cell or a battery. Do not store cells or batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.
- Do not remove a cell or battery from its original packaging until required for use.
- Do not subject cells or batteries to mechanical shock.
- In the event of a cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- Do not use any charger other than that specifically provided for use with the equipment. Refer to the manufacturer's instructions for proper charging instructions.
- Do not use any cell or battery which is not designed for use with the equipment.
- Always purchase the battery recommended by the device manufacturer for the equipment.
- Keep cells and batteries clean and dry.
- Wipe the cell or battery terminals with a clean, dry cloth if they become dirty.
- Do not leave a battery on prolonged charge when not in use.
- After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.
- Retain the original product literature for future reference.
- Use only the cell or battery in the application for which it was intended.
- When possible, remove the battery from the equipment when not in use.

Please remove battery using specified tools. Removing the battery by hand is not allowed. Veuillez retirer la batterie en utilisant les outils

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-
- Dispose of the battery properly.

Disposal

This equipment contains certain materials and chemical compounds incidental to the manufacture of electrical and electronic equipment, and improper “end-of-life” disposal of such equipment can result in environmental contamination. Therefore, this equipment should not be disposed of as ordinary household waste but should instead be delivered to a designated electrical and electronic waste disposal or recycling center. For further information on disposing of electrical and electronic waste, contact the cognizant authority within the local jurisdiction.

Dispose of the Scanner Tips according to standard operating procedures or local regulations for the disposal of contaminated medical waste. For additional Scanner Tips, contact your dealer.

Cleaning, Disinfecting, Sterilizing

Perform the following maintenance activities on your scanner and accessories regularly.

To ensure maximum hygienic safety for the patient, carefully follow the instructions to prepare the scanner for use.

To ensure maximum hygienic safety for the patient and to minimize the risk of cross-contamination, carefully perform the following maintenance activities on your scanner and accessories.

After each patient:

- Clean and disinfect the scanner. See “Clean and Disinfect the Scanner”).
- Clean the Scanner Tip, and then perform Autoclave Sterilization (See “Clean the Scanner Tip” and “Sterilize the Scanner Tip”).

There are 4 models of Scanner Tip:

| Model | Size | UDI-DI | Manual Cleaning | 134°C Sterilization |
|-------|-------|--------------------|-----------------|---------------------|
| TP202 | Large | (01)06973993441034 | Yes | Yes |
| TP203 | Small | (01)06973993441041 | Yes | Yes |
| TP204 | Large | (01)06973993441386 | Yes | Yes |
| TP205 | Small | (01)06973993441393 | Yes | Yes |

Clean and Disinfect the Scanner

General Warnings

WARNING

- Read and follow the warnings and personal protection instructions provided in the Safety Data Sheet (SDS) for the disinfectant used to process the scanner.
- You must wear gloves while cleaning and disinfecting the scanner.
- The scanner must be disinfected with a recommended intermediate-level disinfectant solution with tuberculocidal activity between patients.

-
- DO NOT use a disinfectant containing phenolics or iodophors; doing so will damage the surface coating of the scanner.
 - Never put the scanner in an autoclave device or immerse it in water or the disinfectant solution.
 - Excessive fluids can damage the scanner.
 - Do not use cotton, cloth, or tissues soaked with disinfectant to disinfect the scanner.

Clean the Scanner

- If the scanner is visibly contaminated with blood and/or body fluids, you must clean it before disinfecting it.
- To clean the scanner, follow these steps:
 - ◆ Dampen (do not soak) a lint-free cloth with lukewarm water.
 - ◆ Remove the blood and/or body fluids with the dampened lint-free cloth.

Disinfect the Scanner

After each patient, the scanner must be thoroughly disinfected.

To adequately disinfect the scanner, follow the disinfectant manufacturer's instructions for the appropriate contact time.

CAUTION: If the scanner is visibly soiled, it must be thoroughly cleaned prior to disinfecting. See “Clean the Scanner”

To disinfect the scanner, follow these steps:

- Remove the Scanner Tip.
- Remove all visible soil (see “Clean the Scanner”).
- Use a commercially prepared intermediate level disinfectant wipe. Follow the manufacturer’s instructions for contact time.

Recommended disinfectant wipes: CaviWipes

WARNING: Using a disinfectant that has not been approved may cause damage to the scanner.

- Thoroughly wipe all surfaces of the scanner, DO NOT allow liquid to enter through the gap, air outlet, or pin holes.

WARNING: Do not rinse.

- Allow to air dry.
- After the scanner has dried, use a clean, lint-free cloth dampened with water to remove residual disinfectant from the surface of the scanner.

Clean and Sterilize the Scanner Tip

WARNING

- Wear gloves when handling a contaminated Scanner Tip.
- Read and follow the warnings and personal protection instructions provided in the manufacturer's SDS for the detergent used to clean the Scanner Tip.
- Do not soak the Scanner Tip in disinfectant for a long period.
- Dry the Scanner Tip thoroughly before mounting onto the scanner.
- Do not use an ultrasonic cleaning machine to clean the Scanner Tip.
- Do not soak the Scanner Tip with alcohol-based disinfectants.

Clean the Scanner Tip

To manually cleaning the Scanner Tip, follow these steps:

- Rinse excess soil from the Scanner Tip (2 minutes).
- Using a soft brush, apply enzymatic detergent solution (e.g., Metrex EmPower) to all surfaces.
- Rinse under clean, running water (2 minutes).
- Inspect Scanner Tip. If the Scanner Tip is not clean, repeat the last three steps.
- Use a lens tissue or lint-free cloth to remove any dust from the mirror in the Scanner Tip.
- Dry the scanner tip's mirror carefully with a lens tissue or lint-free cloth.

Sterilize the Scanner Tip

Scanner Tips shipped non-sterile. You must sterilize them before use.

CAUTION: For TP202 and TP203, if you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 60 cycles.

CAUTION: For TP204 and TP205, if you limit the exposure time at 134°C to not more than 6 minutes, you can

autoclave the Scanner Tip up to 180 cycles.

To sterilize the Scanner Tip, follow these steps:

- Put the Scanner Tip into a sealed steam sterilization pouch.
- Place the Scanner Tip in a steam autoclave for sterilization.
 - ◆ Exposure temperature should be set as 134°C.
 - ◆ Exposure time should exceed 3 minutes.
 - ◆ Exposure time should not exceed 6 minutes.

Precautions before use

Perform the following activities on your product and accessories before use.

Visually Inspecting the Scanner for Damage

Visually inspect the scanner for damage or signs of deterioration by doing the following:

- ◆ Inspect the Lens Window.
- ◆ Inspect around the buttons and contact points.

If damage is noted, do not use the scanner and contact your representative.

Visually Inspecting the Scanner Tip

Visually inspect the Scanner Tip for signs of deterioration by doing the following:

- ◆ Verify that the Scanner Tip is not damaged, and its components are not detached.
- ◆ Verify that the Scanner Tip mirror does not have any smudges or scratches on it.

If deterioration is noted, replace the Scanner Tip.

WARNING

- The Lens Window on the scanner is a delicate optical component. Mount the front Protective Cover to protect the Lens Window from damage and dirt when the scanner is not in use.
- The mirror in the Scanner Tip is a delicate optical component. Its clean and undamaged surface is critical to scan quality.

In the event that you see poor scan quality or an unclear video preview in the software, clean the Scanner Tip mirror and the Lens Window of the Handpiece using a microfiber cleaning swab, applying ethanol that is free of impurities.

Troubleshooting

| Problem Description | Action |
|---|--|
| There is mismatching and overlap on the 3D model. | Remove mismatched data and excessive tissue using the Cut tool and rescan. |
| After bite registration, there is a gap or intersection between the upper jaw and the lower jaw. | <ul style="list-style-type: none"> ● Adjust the "Occlusal pressure level" in Preferences, then disable and enable "Occlusal pressure adjustment". ● Delete the incorrect bite view, and rescan. |
| Precision degradation is observed, or images are not well-stitched during acquisition. | Ensure that the Lens Window at the base of the scanner is clean by wiping it with a moist, lint-free cloth or lens tissue. Use a lens tissue or lint-free cloth to remove any dust or water stains. Make sure the Scanner Tip is firmly installed and there are no dark edges on the live video. |
| Reconstruction of metallic preparations is sometimes difficult. | <ul style="list-style-type: none"> ● Adjust the scanner position (for example: distance or angle) and scan more of the area. ● Move the surgical light away from the patient to decrease light scatter. ● Turn on "Shining Surface" button. |
| The Scanner Tip is installed but not detected. No live video is displayed, and the Scanner Tip not detected icon is displayed at the bottom-right of the interface. | Reinstall the Scanner Tip, and make sure the Scanner Tip is in firm contact with the scanner. |
| Fogging appears on the inner surface of the Lens Window at the base of the scanner. | <ul style="list-style-type: none"> ● Mount a completely dry Scanner Tip on the scanner and place the scanner in the holder or set it on the desk, and wait until the fogging fades. If the fogging does not disappear completely after 24 hours, contact your local service provider for assistance. ● Ensure that the Scanner Tip is thoroughly dry before mounting on the scanner, and do not use a cloth soaked in disinfectant to clean the scanner. |

| | |
|---|---|
| <p>The scanner does not cast light, and a static preview image is displayed in the video preview screen.</p> | <p>Close the software, and then launch the software again.</p> |
| <p>Scanner is not shown in the 'Scanner List' in the software while the Status Indicator keeps blinking in breathing mode</p> | <p>Reset the network in the scanner. To do so, while the Status Indicator is in breathing mode, simultaneously press the Scan Button on both sides for two seconds and release. The Status Indicator should start blinking rapidly, then enter breathing mode. The device should then be shown in the 'Scanner List' in the software.</p> |

Regulatory Information

Regulatory Information

The device complies with the following regulations:

| |
|--|
| MDR: (EU) 2017/745 Medical Device Regulation, Class I following the Rule 5. |
| FDA Center for Devices & Radiological Health CDRH - Title 21 CFR 872.3661 (USA). |
| Medical Devices Regulations (Canada). |
| RoHS: Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment, 2011/65/EU Annex II and its amendment Directive (EU) 2015/863 |
| RED: Directive 2014/53/EU The Radio Equipment Directive |
| FCC: Part 15 of The Federal Communications Commission Rules |
| ISED: Innovation, Science and Economic Development Canada |

Compliance with European and International Standards

EN / IEC 60601-1: Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance

ANSI/AAMI ES 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

CAN/CSA-C22.2 No. 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

EN / IEC 60601-1-2: Medical Electrical Equipment, Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests

EN / IEC 80601-2-60: Medical electrical equipment — Part 2-60: Particular requirements for the basic safety and essential performance of dental equipment

EN / IEC 62471: Photobiological safety of lamps and lamp systems

EN / ISO 17664: Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices

EN / ISO 17665-1: Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

EN / IEC 60601-1-6: Medical Electrical Equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability

EN / IEC 62366-1: Medical devices - Part 1: Application of usability engineering to medical devices

EN / IEC 62304: Medical device software - Software life cycle Processes

EN ISO 10993: Biological evaluation of medical devices

ISO 14971: Medical devices - Application of risk management to medical devices

EN / ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements

EN ISO 20417: Medical devices — Information to be supplied by the manufacturer

ISO 9687: Dentistry - Graphical symbols for dental equipment

AAMI TIR 12: Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers

AAMI TIR 30: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices

EN / IEC 62133-2: Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems

EN 50566: Product standard to demonstrate the compliance of wireless communication devices with the basic restrictions and exposure limit values related to human exposure to electromagnetic fields in the frequency range from 30 MHz to 6 GHz: hand-held and body mounted devices in close proximity to the human body.

EN 301489-1: Electromagnetic compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements; Harmonized Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU and the essential requirements of article 6 of Directive 2014/30/EU

EN 301489-17: Electromagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems; Harmonized Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU

EN 301893: 5GHz RLAN; Harmonized Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU

Classification in Accordance with EN/IEC 60601-1

Type of protection against electric shock: Internally powered

Degree of protection against electric shock: Type BF Applied Part

Mode of operation: Continuous operation

Flammable anesthetics: Not suitable for use in the presence of flammable anesthetics or a mixture of flammable anesthetics with air or oxygen or nitrous oxide.

Conformity with EN/IEC 60601-1-2

IEC 60601-1-2 EMC requirements and tests, Medical Electrical Equipment including CISPR 11 Group 1, Class B.

Electromagnetic Compatibility Precautions

Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC). Medical equipment must be installed and put into service according to the EMC information provided in this documentation.

Other equipment can interfere with communications with the device, even if the equipment complies with CISPR emissions requirements.

WARNING: Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the DEVICE, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could happen.

WiFi

The device operates with the 802.11a/n/ac protocol. Channel 38 or 46 is priority used for the handpiece. The channel bandwidth is 40MHz. The frequency range is 5150-5250MHz,5725-5850MHz (the actual frequencies are dependent on local regulations and the configuration of the product). The maximum output power is 17.88 dBm.

This device complies with part 15 of the FCC Rules and contains license exempt transmitter(s)/receiver(s) that comply with ISSED's license-exempt RSS(s).

Operation is subject to the following two conditions:

- This device may not cause interference.
- This device must accept any interference, including interference that may cause undesired operation of the device.

NOTE: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

SAR Value for Handpiece:

0.798 W/Kg,10g for CE

0.011 W/Kg,1g for FCC

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Limited by local law regulations, the version for North America does not have a region selection option.

The device is for indoor use only and operates in the 5150-5250 MHz band to reduce the potential for harmful interference to co-channel mobile satellite systems.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not in-stalled and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Restrictions in the 5 GHz band: According to Article 10 (10) of Directive 2014/53/EU, the packaging shows that this radio equipment will be subject to some restrictions when placed on the market in Belgium (BE), Bulgaria (BG), the Czech Republic (CZ), Denmark (DK), Germany (DE), Estonia (EE), Ireland (IE), Greece (EL), Spain (ES), France (FR), Croatia (HR), Italy (IT), Cyprus (CY), Latvia (LV), Lithuania (LT), Luxembourg (LU), Hungary (HU), Malta (MT), Netherlands (NL), Austria (AT), Poland (PL), Portugal (PT), Romania (RO), Slovenia (SI), Slovakia (SK), Finland (FI), Sweden (SE), Turkey (TR), Norway (NO), Switzerland (CH), Iceland (IS), and Liechtenstein (LI). The WLAN function for this device is restricted to indoor use only when operating in the 5150 to 5250 MHz frequency range.

| | |
|---|--|
|  | AT, BE, BG, CH, CY, CZ, DE, DK, EE, EL, ES, FI, FR, HR, HU, IE, IS, IT, LI, LT, LU, LV, MT, NL, NO, PL, PT, RO, SE, SI, SK, TR |
|---|--|

Guidance and Manufacturer’s Declaration

Guidance and Manufacturer’s Declaration - Electromagnetic Emission

The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

| Emissions Test | Compliance | Electromagnetic Environment - |
|----------------|------------|-------------------------------|
|----------------|------------|-------------------------------|

| | | Guidance |
|--------------------------|--------------------|---|
| RF Emissions CISPR 11 | Group 1 Class B | The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems (from Non-RF wireless communication equipment)

The DEVICE is intended for use in the electromagnetic environment specified below (IEC 60601-1-2 Table 4). The customer or the user of the DEVICE should assure that it is used in such an environment.

| Immunity Test | IEC 60601 Test Level | Compliance Level | Electromagnetic Environment - Guidance |
|--|-----------------------------|-----------------------------|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±15 kV air | ±8 kV contact ±15 kV air | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| Radiated RF IEC 61000-4-3 | 3 V/m 80MHz – 2.7GHz | 3 V/m 80MHz – 2.7GHz | Environment of a professional healthcare facility. |

NOTE: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DEVICE is used exceeds the applicable RF compliance level above, the DEVICE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the DEVICE.

Guidance and Manufacturer’s Declaration – Electromagnetic Immunity for Equipment and Systems (from RF wireless communication equipment)

For the immunity to proximity fields from RF wireless communications equipment, the DEVICE is compliant with the test levels specified below, according to IEC60601-1-2 Table 9. The customer or user of the DEVICE should assure that it is used in such an environment.

| Test Frequency (MHz) | Band (MHz) | Immunity Test Levels |
|----------------------|-------------|--|
| 385 | 380 – 390 | Pulse modulation 18Hz, 27V/m |
| 450 | 430 – 470 | FM, ± 5 kHz deviation, 1 kHz sine, 28V/m |
| 710 | 704 – 787 | Pulse modulation 217Hz, 9V/m |
| 745 | | |
| 780 | | |
| 810 | 800 – 960 | Pulse modulation 18Hz, 28V/m |
| 870 | | |
| 930 | | |
| 1720 | 1700 – 1990 | Pulse modulation 217Hz, 28V/m |
| 1845 | | |
| 1970 | | |
| 2450 | 2400 – 2570 | Pulse modulation 217Hz, 28V/m |
| 5240 | 5100 – 5800 | Pulse modulation 217Hz, 9V/m |
| 5500 | | |
| 5785 | | |

Accessories

The use of accessories other than those specified, except for those sold by the manufacturer of the equipment, as

replacement parts for internal components may result in increased emissions or decreased immunity of the medical equipment.

Other Equipment

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify normal operation.

Technical Specifications

Model

Owandy-IOS W

Technical Specifications

| Components | Technical Specifications |
|------------------------|---|
| Weight | Scanner (including battery, without Scanner Tip): 245 g Tip Large - Sterilization (TP202): 14g Tip Small - Sterilization (TP203): 12g Tip Large - Sterilization (TP204): 14g Tip Small - Sterilization (TP205): 12g |
| Color | 3D full color |
| USB Wireless Card Port | USB 3.0 |
| USB Charge Port | Type C Port |
| Field of View | Tip Large - Sterilization (TP202): 16mm x 14mm Tip Small - Sterilization (TP203): 12mm x 12mm Tip Large - Sterilization (TP204): 16mm x 14mm Tip Small - Sterilization (TP205): 12mm x 12mm |

| | |
|---|---|
| Minimum System Requirements of Laptop | Processor: 10th Generation Intel® Core™ i7-10750H or AMD Ryzen™ 7 4700U Memory: 16 GB Graphics card: NVIDIA GeForce GTX 1650 4GB Disk: 512 GB SSD Operating system: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 port Optional: Touch screen |
| Minimum System Requirements of Desktop | Processor: 9th Generation Intel® Core™ i5-9600 or AMD Ryzen™ 5 2600 Memory: 16 GB Graphics card: NVIDIA GeForce GTX 1650 4GB Disk: 512 GB SSD Operating system: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 port Optional: Touch screen |

CAUTION: It is MANDATORY to check that your system configuration is compatible with the computer system requirements for the accompanying software.

Environmental Requirements

| Components | Environmental Requirements |
|---|----------------------------|
| Operating Temperature | 15°C ~ 30°C |
| Transport and Storage Temperature | -10°C ~ 50°C |
| Operating Relative Humidity | 10% ~ 65% RH |
| Transportation and Storage Relative Humidity | 10% ~ 95% RH |
| Operating Atmospheric Pressure | 70 ~ 106 KPa |
| Transportation and Storage Atmospheric Pressure | 60 ~ 106 KPa |

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